



DERMA GRIP HIDE PASTE

#1232 Taxidermy Adhesive and Hide Paste is suitable for all tanning methods and all animal skins. Allows up to six hours of working time. Forms a strong, final bond with excellent holding strength and slip resistance. Ideal for dry-preserved tans.

91-1232A	12 oz spout pack	\$ 22.00
91-1232G	1 Gallon	\$ 65.00
91-1232G2	2 Gallons	\$ 125.00



DERMA GRIP DRY HIDE PASTE

#1232 Dry Taxidermy Adhesive and Hide Paste is suitable for all tanning methods and all animal skins. Forms a strong, final bond with excellent holding strength and slip resistance. Mixing only amounts as needed allows for a longer shelf life. Ideal for dry-preserved tans. Mix with water to a 1:1 ratio.

91-12325	5 lbs	\$ 65.00
----------	-------	----------



JAW JUICE

97-JJ100	3.7 Ozs	\$ 13.25
----------	---------	----------



HIDE GLUE

Save time with this pre-mixed thin spread adhesive. Apply to manikin with a paint brush. Allows adequate working time. No carding necessary, just work out any air pockets the first couple of days. Easily cleans up with water.

91-GL10	1 Gallon	\$ 50.00
---------	----------	----------



RON SCHAEFER HIDE GLUE

This is a water based clay glue. This glue has been specially formulated for the taxidermy industry. The glue has great adhesion with the flexibility to move and adjust the skin for days. This re-hydrate able glue is specially designed for African animals but works incredible for North American animals as well. Freezable/thaw able.

91-HG10	1 Gallon	\$ 66.00
91-HG50	5 Gallons (DS)	\$ 225.00



EPO-GRIP WATERBASE HIDE ADHESIVE

Use this high performance product where extended working times are required with little pinning or carding. Pre-colored pink tone. Easily applied by glove or brush. Clean up easily with a damp cloth. INTENDED FOR TANNED HIDE APPLICATION ONLY.

91-EGWHQ	Quart	\$ 21.75
91-EGWHG	Gallon (DS)	\$ 56.00



HIDE PASTE

Allows you to move and adjust your skin on the manikin. Once dry it adheres the skin to the manikin, retaining the details. Mix with warm water to form a paste and apply to manikin.

90-HP05	5 lbs	\$ 16.00
90-HP10	10 lbs	\$ 32.00
90-HP25	25 lbs	\$ 81.00



DEXTRIN

Use to make various pastes and machés. Add to plaster for a binder to make it harden.

90-DX05	5 lbs	\$ 31.00
90-DX10	10 lbs	\$ 62.00



THE ORIGINAL B.S. HIDE PASTE (BRYAN'S SUPREME)

The industry's leading Hide Paste for over 25 years. Now reformulated so that it also works with dry preservative.

91-BS1Q	Quart	\$ 20.00
91-BS1G	1 Gallon	\$ 57.00



HIDE PASTES / ADHESIVES



EPO-GRIP HIDE ADHESIVE (EPOXY)

This high strength system produces superior results. When applied to a form it does not sag, slip, run or shrink. Cures in 4-6 hours. Saves time, pinning is not usually required. Works for all tanning methods and all animals.

91-EGHAQ Quart \$ 38.00
 91-EGHAG Gallon \$ 104.50



EPO-GRIP #30 PASTE (Slow Set)

An opalescent, two component, light paste, epoxy resin adhesive. High strength and minimum shrinkage. Initial cure occurs in 3-4 hours. Does not run or sag. Use as hide adhesive for fish and fish seam sealer or anywhere a slow set high strength adhesive is needed.

91-EG3008 8 ozs \$ 30.50
 91-EG3016 16 ozs \$ 44.00



EPO GRIP CLEAR PASTE

A clear, two component, fast setting, epoxy paste adhesive which is easily mixed and initial cure occurs in 6-10 minutes. No shrinkage. Can be cut in 15 minutes and sanded in 30 minutes. Use to install ear liners, fish fins, artificial heads, fix fins, fin backer, seam sealer and general shop repairs.

91-EGCL08 8 ozs \$ 30.50
 91-EGCL16 16 ozs \$ 44.00



EPO-GRIP EAR MAGIC

A time saving reinforced translucent ear adhesive designed to produce strong, flexible ear liner quality ears without the mess. Can also be used to install ear liners.

91-EGEM08 8 ozs \$ 33.25
 91-EGEM16 16 ozs \$ 46.75
 91-EGEM32 32 ozs (DS) \$ 77.00



EPO-GRIP MACHE MAGIC

An innovative time saving filler. This two component system is easily mixed and initial curing occurs in 6-10 minutes (77° F). Mache Magic can be easily tinted. This material can be used regardless if the piece is wet or dry. With latex gloves and the aide of a little soap and water the material can be smoothed and feathered to reducing finish work in transition areas. Works great for replacing cheek flesh, building gum lines, tongues, installing jaw sets and rebuilding transition areas.

91-EGMM08 8 ozs (DS) \$ 26.50
 91-EGMM16 16 ozs (DS) \$ 39.00
 91-EGMM32 32 ozs (DS) \$ 68.25



EPO-GRIP PASTE FILLER

A specially formulated two component paste consisting of 100% epoxy solids and bonding agents. Being neutral in color, this filler can be easily tinted. Epo-Grip Paste Filler is easily mixed and initial curing occurs in 6-10 minutes (at 77 degrees F). This filler sands, drills and taps within a 1 hour period. Use to repair skulls, teeth, bone antlers & horns, filling manikin repairs and embedding bolts, screws, etc.

91-EGPA08 8 ozs (DS) \$ 26.50
 91-EGPA16 16 ozs (DS) \$ 39.00
 91-EGPA32 32 ozs (DS) \$ 68.25



ELIMINATES
 SANDING OR
 SCRATCHING
 OF FORMS.

MANNIKIN PREP

Mannikin Prep is specially formulated for use as an adhesive primer coat between any hide paste and all form surfaces. Can be used on all forms, earliners and as an advanced adhesive for habitat and base work.

91-MP16 16.75 ozs \$ 21.00



MAGIC SMOOTH

This two-part epoxy has a creamy consistency. Use it to fill in seams or shrinkage, smooth out rough areas, repair fins, hide paste or installing ear liners. Semi-transparent.

91-WMSP Pint \$ 32.00
 91-WMSQ Quart \$ 82.50



EAR ADHESIVE FAST

Fast setting two-part epoxy. Mix equal amounts resin and hardener. Spread a thin layer on the ear liner and the skin. Apply pressure. Cure time is 10-20 minutes. May also be used as an all-purpose adhesive. Can be tinted. Dries semi-transparent.

91-WEARP Pint \$ 32.00



Zap Products



ZAP A GAP

This medium viscosity instant glue allows approximately 7-10 seconds working time with an average of 20 seconds to cure. Can use Zap Kicker to accelerate cure time. Will bond virtually anything to anything. Even works on oily surfaces. Superior gap filling formula.

91-ZAP1	1 oz	\$ 9.25
91-ZAP5	(5) .5 grams	\$ 6.00



ZAP KICKER

The best way to accelerate the cure time of all brands of super glues. When sprayed directly on a glued joint it forces the immediate cure of the super glue.

91-ZKICK	2 ozs	\$ 8.00
----------	-------	---------



ZAP GEL

Zap Gel is the thickest and strongest super glue available. Will not run or flow. Can be applied to vertical or overhead applications. The tube comes with a long tapered tip allowing a precise application. Use Zip Kicker to speed the cure time.

91-ZGEL03	3 Grams	\$ 2.00
91-ZGEL20	1 oz	\$ 11.00



CURVED TIP SYRINGE

92-HPS412	12 cc	\$ 1.50
-----------	-------	---------



FLEXI-TIPS

Will fit most brands of super glue bottles. Allows for smaller drops. Extends into hard to reach places with precise application.

91-ZTIP	12 Tips	\$ 2.75
---------	---------	---------



HYPODERMIC SYRINGES

92-HPS12	12 cc	\$ 1.25
92-HPS20	20 cc	\$ 2.00
92-HPS35	35 cc	\$ 2.50
One (1) needle included with each syringe.		
92-HPSNED	Extra Needle	\$.50



ADVANCED TANNING SOLUTIONS INJECTION FLUID

For injecting bird feet, small mammal toes and antlers in velvet to preserve and eliminate shrinkage.

91-AJECTQ	Quart	\$ 28.00
91-AJECTG	Gallon (DS)	\$ 76.50



FANTASTIC CAST INJECTION GRADE

A two part liquid used to inject bird feet or foot pads of small mammals to prevent shrinkage and help keep their shape. Can also be used to fill in drumming. Sets in 12-18 hours.

91-FCIG4	(2) 4 ozs Each	\$ 24.00
----------	----------------	----------



BIRD FEET INJECTION FLUID

A preservative for injection into bird feet, wings, feet of small game and velvet antlers to reduce shrinkage. Safe to use without offensive odors.

91-K130P	Pint	\$ 23.00
91-K130Q	Quart	\$ 37.25



EPO GRIP LIQUID FAST SET

A clear, low viscosity, liquid epoxy adhesive. Initial bonding occurs in 6-10 minutes. Low shrinkage. Use to install ear liners, fish fins, artificial heads, eyes and general shop repairs.

91-EGLQ08	8 ozs	\$ 39.75
91-EGLQ16	16 ozs (DS)	\$ 55.50



Injectables